


## Anti-Mono-methyl-Histone H4 (Arg3) Polyclonal Antibody 2

### Product Details

Ig Type: IgG  
Reactivity: Human  
Conjugation: Unconjugated  
Molecular Weight: Theoretical: 14 kDa.

### Applications

Verified Activity: 1. Western blot analysis of extracts from Hela cells, 1:2000.



Application: WB  
Recommended WB: 1:1000-2000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: A synthesized monomethylpeptide: Histone H4 around the monomethylation site of Arg3  
Synonyms: MonoMe-H4R3; Mono-Me-Histone H4 (Arg3); Mono-methyl-Histone H4 (R3); H4R3me; Mono-Me-Histone H4 (R3); Histone H4R3-monomethylated

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481