


Anti-Mono-methyl-Histone H4 (Lys79) Polyclonal Antibody 2

Product Details

Ig Type: IgG
Reactivity: Human,Mouse,Rat
Conjugation: Unconjugated
Molecular Weight: Theoretical: 14 kDa.

Applications

Verified Activity: 1. Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:2000.



Application: WB
Recommended WB: 1:500-1000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized monomethylpeptide: Histone H4 around the monomethylation site of Lys79
Synonyms: Mono-methyl-Histone H4 (K79);Mono-Me-Histone H4 (Lys79);H4K79me;Histone H4K79-monomethylated;Mono-Me-Histone H4 (K79);MonoMe-H4K79

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481