


## Anti-Phospho-Histone H4 (Ser47) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 14 kDa.
Purification:	Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

## Applications

Verified Activity: 1. Western blot analysis of extracts from Hela cells, untreated (-) or treated with Nocodazole (1µg/ml, 24 hr; +), using TMAC-01969 diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:1000-2000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Peptide sequence around phosphorylation site of serine 47(R-I-S(p)-G-L) derived from Human Histone H4
Antigen Species:	Human
Uniprot ID:	P62805
Synonyms:	Histone H4 (p-Ser47);Histone H4 (p-S47);p-Histone H4 (Ser47);p-Histone H4 (S47)

## Research Background

Histone H4 is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells. H4 is a structural component of the nucleosome, and is subject to covalent modification, including acetylation and methylation, which may alter expression of genes located on DNA associated with its parent histone octamer.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481