

Anti-Phospho-HNF1 α (Ser247) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Rat,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 69 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Applications

Application:	WB
--------------	----

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human HNF1 α around the phosphorylation site of Ser247
Antigen Species:	Human
Gene ID:	6927
Uniprot ID:	P20823
Synonyms:	HNF1 α (p-S247);p-HNF1 α (Ser247);p-HNF1 α (S247);HNF1 α (p-Ser247)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481