

Mouse Anti-Human IgG4 Secondary Antibody-HRP (6Z129)

Product Details

Conjugation:	Unconjugated
Clone:	6Z129
Purification:	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.

Applications

Application:	ELISA,IF,IHC,WB
Recommended	WB: 1:1000-10000 ELISA 1:10000-20000 IHC: 1:50-100 IF: 1:50-1000

Properties

Stability & Storage:	Store at -15°C to -25°C for 12 months (Do not lower than -25°C). Avoid repeated freeze-thaw cycles. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481