

Anti-Hsp27/HSPB1 Antibody (9E653)



Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 27 kDa.
Clone:	9E653

Applications

Verified Activity:

1. Western blot analysis of Hela with HSP27 Mouse mAb diluted at 1:2,000.
2. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using HSP27 Mouse mAb diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:1000-3000; IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Uniprot ID: P04792

Synonyms: HSP28;HMN2B;CMT2F;Hsp27;HSPB1;HSP 27;SRP27

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481