

Anti-Human lambda light chain Antibody (8V560)

Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 25 kDa.
Clone:	8V560
Purification:	Affinity purified

Applications

Verified Activity: 1. Western blot analysis on human IgG using anti-human lambda light chain antibody.



Application:	WB
Recommended	WB: 1:1000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	Human IgG
Antigen Species:	human

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481