


## Anti-Phospho-IL-13RA1 (Tyr405) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 49 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western Blot analysis of 3T3 cells using Phospho-IL-13R $\alpha$ 1 (Y405) Polyclonal Antibody.



Application:	ELISA,IF,WB
Recommended	WB: 1:500-2000; IF: 1:200-1000; ELISA: 1:10000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized phosphopeptide: human IL-13R $\alpha$ 1 around the phosphorylation site of Y405
Antigen Species:	Human
Uniprot ID:	P78552
Synonyms:	p-IL-13RA1 (Y405);IL-13RA1 (p-Y405);p-IL-13RA1 (Tyr405);IL-13RA1 (p-Tyr405)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481