

## Anti-Phospho-IRAK4 (Thr345, Ser346) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 60 kDa.
Purification:	Affinity purification using immunogen.

### Applications

Application:	WB
Recommended	WB: 1:1000-3000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human IRAK4 around the phosphorylation site of Thr345/Ser346
Antigen Species:	Human
Uniprot ID:	Q9NWZ3
Synonyms:	IRAK4 (p-T345, S346);p-IRAK4 (Thr345, Ser346);p-IRAK4 (T345, S346);IRAK4 (p-Thr345, Ser346)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481