

Anti-Phospho-IRAK4 (Thr345, Ser346) Polyclonal Antibody

Product Details

| | |
|-------------------|--|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 60 kDa. |
| Purification: | Affinity purification using immunogen. |

Applications

| | |
|--------------|-----------------|
| Application: | WB |
| Recommended | WB: 1:1000-3000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|---|
| Immunogen: | A synthesized phosphopeptide: human IRAK4 around the phosphorylation site of Thr345/Ser346 |
| Antigen Species: | Human |
| Uniprot ID: | Q9NWZ3 |
| Synonyms: | IRAK4 (p-T345, S346);p-IRAK4 (Thr345, Ser346);p-IRAK4 (T345, S346);IRAK4 (p-Thr345, Ser346) |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481