

Anti-Phospho-IRS-1 (Ser794) Polyclonal Antibody



Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 180 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity:

1. Western Blot analysis of Rat cells using Phospho-IRS-1 (S794) Polyclonal Antibody.
2. Western Blot analysis of various cells using Phospho-IRS-1 (S794) Polyclonal Antibody.



Application:	ELISA,IF,IHC,WB
Recommended	WB: 1:500-2000; IHC: 1:100-300; IF: 1:200-1000; ELISA: 1:40000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human IRS-1 around the phosphorylation site of S794
Antigen Species: Human
Uniprot ID: P35568
Synonyms: p-IRS1 (S794);p-IRS1 (Ser794);IRS-1;Insulin receptor substrate 1;IRS1 (p-S794);IRS1 (p-Ser794)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481