


Anti-Phospho-IRS-1 (Tyr896) Polyclonal Antibody

Product Details

| | |
|-------------------|---|
| Ig Type: | IgG |
| Reactivity: | Human,Mouse,Rat |
| Conjugation: | Unconjugated |
| Molecular Weight: | Actual: 180 kDa. |
| Purification: | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |

Applications

Verified Activity: 1. Western blot analysis of KB MCF-7 using p-IRS-1 (Y896) antibody.



Application: ELISA,WB
Recommended WB: 1:500-2000; ELISA: 1:40000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human IRS-1 around the phosphorylation site of Y896
Antigen Species: Human
Uniprot ID: P35568
Synonyms: IRS-1;IRS1 (p-Y896);p-IRS1 (Y896);IRS1 (p-Tyr896);Insulin receptor substrate 1;p-IRS1 (Tyr896)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481