

Anti-Phospho-Mer/TYRO3 (Tyr753/685) Polyclonal Antibody


Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 110 kDa.
Purification:	The antibody was purified by immunogen affinity chromatography.

Applications

Verified Activity:

1. Western blot analysis of MERTK (phospho-Tyr753/685) expression in HEK293T, H1688 whole cell lysates.



Application:	WB
Recommended	WB: 1:500-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human MER/TYRO3 around the phosphorylation site of Tyr753/Tyr685
Antigen Species:	Human
Synonyms:	Mer/TYRO3 (p-Tyr753/685);p-Mer/TYRO3 (Tyr753/685);Mer/TYRO3 (p-Y753/685);p-Mer/TYRO3 (Y753/685)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481