


Anti-Phospho-HGFR/c-Met (Tyr1234, 1235) Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: Human,Mouse,Rat
Conjugation: Unconjugated

Applications

Verified Activity: 1. Western blot analysis of extracts from HUVEC cells untreated or treated with EGF using Met (Phospho-Tyr1234/1235) Antibody TMAC-02566.



Application: WB

Properties

Shipping: Shipping with blue ice.

Antigen Details

Synonyms: HGFR/c-Met (p-Tyr1234, 1235);p-HGFR/c-Met (Tyr1234, 1235);p-HGFR/c-Met (Y1234, 1235);
HGFR/c-Met (p-Y1234, 1235)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481