

Anti-MMP-9 Antibody (7Q238)




Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 78 kDa.
Clone:	7Q238
Purification:	Affinity purification

Applications

Verified Activity:

1. Western Blot analysis using MMP-9 Monoclonal Antibody against HEK293 (1), MMP7-hlgGfc transfected HEK293 (2) cell lysate and MMP9-hlgGfc transfected HEK293 (3) cell lysate.
2. Immunohistochemistry analysis of paraffin-embedded brain tissues with DAB staining using MMP-9 Monoclonal Antibody.
3. Immunofluorescence analysis of NIH/3T3 cells using MMP-9 Monoclonal Antibody (green).



Application: ELISA,FCMM,IF,IHC,WB

Recommended WB: 1:500-2000; IHC: 1:200-1000; IF: 1:200-1000; FCM: 1:200-400; ELISA: 1:10000;

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein: human MMP-9
Antigen Species: human
Gene ID: 4318
Uniprot ID: P14780
Synonyms: GELB;MANDP2;CLG4B;MMP-9;Gelatinase B

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481