

Anti-mStrawberry Polyclonal Antibody

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of mStrawberry recombinant protein, using TMAC-02668 diluted at 1:5,000.



Application:	IF,IHC,IP,WB
Recommended	WB: 1:2000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Synonyms: mStrawberry

Research Background

mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481