


Anti-MUC2 (ABT-MUC2) Antibody (60821)

Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Clone:	60821

Applications

Verified Activity: 1. Immunohistochemical analysis of paraffin-embedded Colon.



Application:	IF,IHC-P
Recommended	IHC-P: 1:100-500

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human MUC2
Antigen Species:	Human
Synonyms:	MUC 2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481