


Anti-MYC Tag Antibody (4B356)

Product Details

| | |
|---------------|--|
| Conjugation: | Unconjugated |
| Clone: | 4B356 |
| Purification: | The antibody was affinity-purified from antiserum by affinity-chromatography |

Applications

Verified Activity: 1. Western blot analysis of His-MBP-C-Myc recombinant protein, diluted at 1) 1:5000 2) 1:10000.



| | |
|--------------|------------------|
| Application: | WB |
| Recommended | WB: 1:5000-10000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|-----------|---------|
| Synonyms: | BHLHE39 |
|-----------|---------|

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481