

## Anti-MYC Tag Antibody (85498)

### Product Details

|              |              |
|--------------|--------------|
| Ig Type:     | IgG1         |
| Conjugation: | Unconjugated |
| Clone:       | 85498        |

### Applications

Verified Activity: 1. Western blot analysis of c-Myc Recombinant protein, diluted at 1:2000.



|              |                |
|--------------|----------------|
| Application: | WB             |
| Recommended  | WB: 1:500-2000 |

### Properties

|                      |   |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping:            | Shipping with blue ice.   |

### Antigen Details

|             |                       |
|-------------|-----------------------|
| Immunogen:  | A synthesized peptide |
| Uniprot ID: | P01106                |
| Synonyms:   | bHLHe39               |

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481