


Anti-MYC Tag Antibody (3I281)

Product Details

Ig Type:	IgG1
Conjugation:	Unconjugated
Clone:	3I281

Applications

Verified Activity: 1. Western blot analysis of His-MBP-C-Myc recombinant protein, diluted at 1) 1:5000 2) 1:10000.



Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide
Uniprot ID:	P01106
Synonyms:	bHLHe39

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481