

Anti-MPO Antibody (8Z615)



Product Details

Ig Type:	IgG2a, Kappa
Reactivity:	Human,,Mouse,Rat
Conjugation:	Unconjugated
Clone:	8Z615
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Applications

Verified Activity:

1. HL60 whole cell lysates were separated by 12% SDS-PAGE, and the membrane was blotted with anti-MPO (ABT-MPO) antibody.
2. Human spleen tissue was stained with Anti-Myeloperoxidase (ABT-MPO) Antibody.



Application:	IF,IHC,WB
Recommended	WB: 1:200-1000; IHC: 1:100-500; IF: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide: human Myeloperoxidase(MPO)
Antigen Species: Human
Gene ID: 4353
Uniprot ID: P05164
Synonyms: myeloperoxidase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481