

## Anti-Nano Tag9 Antibody (1T795)

### Product Details

Ig Type: IgG1  
Conjugation: Unconjugated  
Clone: 1T795

### Applications

Verified Activity: 1. Western blot analysis of Recombinant Nano-Tag9 Protein with Nano-Tag9 Mouse mAb diluted at 1)1:10,000 2)1:5,000.

4CL



Application: WB  
Recommended WB: 1:2000-5000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: A synthesized peptide  
Synonyms: Nano Tag 9

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481