

## Anti-Phospho-NBK (Ser35) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 18 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

### Applications

Application:	ELISA,IF,IHC-P
Recommended	IHC-P: 1:100-300; ELISA: 1:5000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human NBK around the phosphorylation site of Ser35
Antigen Species:	human
Uniprot ID:	Q13323
Synonyms:	NBK (p-S35);p-NBK (Ser35);p-NBK (S35);NBK (p-Ser35)

### Research Background

**BCL2 interacting killer(BIK)** Homo sapiens The protein encoded by this gene shares a critical BH3 domain with other death-promoting proteins, such as BID, BAK, BAD and BAX, that is required for its pro-apoptotic activity, and for interaction with apoptotic members of the BCL2 family, and viral survival-promoting proteins. Since the activity of this protein is suppressed in the presence of survival-promoting proteins, it is suggested as a likely target for apoptotic proteins. [provided by RefSeq, Sep 2011],

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481