

Anti-Acetyl-NF-E4 (Lys43) Polyclonal Antibody



Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 19 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity:

1. Western Blot analysis of AD-293 cells using Acetyl-NF-E4 (K43) Polyclonal Antibody.
2. Western Blot analysis of AD-293 cells using Acetyl-NF-E4 (K43) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:20000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized acetylpeptide: human NF-E4 around the acetylation site of K43
Antigen Species: Human
Uniprot ID: P60484
Synonyms: Transcription factor NF-E4;Acetyl-NF-E4 (K43);Ac-NF-E4 (K43);NF-E4K43ac;Ac-NF-E4K43;Ac-NF-E4 (Lys43);NF-E4K43-acetylated;NFE4

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481