

Anti-NFKB1 Antibody (3N127)




Product Details


Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 50/120 kDa.
Clone:	3N127
Purification:	Affinity-chromatography

Applications

Verified Activity:

1. Western blot analysis of NF- κ B p105/p50 expression in (1)PC-12 cell lysate; (2) NIH/3T3 cell lysate.
2. All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.
3. All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.
4. All lanes use the Antibody at 1:3K dilution for 1 hour at room temperature.
5. Immunohistochemical analysis of paraffin-embedded human kidney, using NF- κ B p105/p50 Antibody.





Application: IHC, WB

Recommended WB: 1:500-2000; IHC: 1:50-200

Properties

Stability & Storage: Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide: human NF- κ B (p105/p50)

Antigen Species: Human

Uniprot ID: P19838

Synonyms: Nuclear factor NF-kappa-B p105 subunit; DNA-binding factor KBF1; EBP-1; NFKB1; Nuclear factor of kappa light polypeptide gene enhancer in B-cells 1

Research Background

NF- κ B-p105 a transcription factor of the nuclear factor-kappaB (NF- κ B) group. Undergoes cotranslational processing by the 26S proteasome to produce a 50 kD protein. The 105 kD protein is a Rel protein-specific transcription inhibitor and the 50 kD protein is a DNA binding subunit of NF- κ B. NF- κ B is a transcription regulator that is activated by various intra- and extra-cellular stimuli such as cytokines, oxidant-free radicals, ultraviolet irradiation, and bacterial or viral products.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481
