




## Anti-Phospho-NF-kB p65 (Thr435) Polyclonal Antibody 2

### Product Details

Ig Type:	IgG
Reactivity:	Human, Mouse, Rat
Conjugation:	Unconjugated
Purification:	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

### Applications

- Verified Activity:
- Western blot analysis of extracts from HeLa cells untreated or treated with TNF $\alpha$ +CA using NFkB-p65 (Phospho-Thr435) Antibody TMAC-02816.
  - Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using NF-kB p65 (phospho-Thr435) antibody (TMAC-02816).
  - Immunofluorescence staining of methanol-fixed HeLa cells using NF-kB p65 (phospho-Thr435) antibody (TMAC-02816, Red).



Application: IF, IHC, WB

### Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

---

### Antigen Details

Immunogen:	Peptide sequence around phosphorylation site of threonine 435(E-G-T(p)-L-S) derived from Human NFκB-p65
Antigen Species:	Human
Uniprot ID:	Q04206
Synonyms:	NFKB3;p65;NF-κB p65 (p-T435);p-NF-κB p65 (Thr435);v-rel avian reticuloendotheliosis viral oncogene homolog A;p-NF-κB p65 (T435);NF-κB p65 (p-Thr435)

---

### Research Background

NF-kappa-B is a pleiotropic transcription factor which is present in almost all cell types and is involved in many biological processes such as inflammation, immunity, differentiation, cell growth, tumorigenesis and apoptosis. NF-kappa-B is a homo- or heterodimeric complex formed by the Rel-like domain-containing proteins RELA/p65, RELB, NFKB1/p105, NFKB1/p50, REL and NFKB2/p52 and the heterodimeric p65-p50 complex appears to be most abundant one. The dimers bind at kappa-B sites in the DNA of their target genes and the individual dimers have distinct preferences for different kappa-B sites that they can bind with distinguishable affinity and specificity. Different dimer combinations act as transcriptional activators or repressors, respectively. NF-kappa-B is controlled by various mechanisms of post-translational modification and subcellular compartmentalization as well as by interactions with other cofactors or corepressors. NF-kappa-B complexes are held in the cytoplasm in an inactive state complexed with members of the NF-kappa-B inhibitor (I-kappa-B) family. In a conventional activation pathway, I-kappa-B is phosphorylated by I-kappa-B kinases (IKKs) in response to different activators, subsequently degraded thus liberating the active NF-kappa-B complex which translocates to the nucleus. NF-kappa-B heterodimeric p65-p50 and p65-c-Rel complexes are transcriptional activators. The NF-kappa-B p65-p65 complex appears to be involved in invasion-mediated activation of IL-8 expression. The inhibitory effect of I-kappa-B upon NF-kappa-B in the cytoplasm is exerted primarily through the interaction with p65. p65 shows a weak DNA-binding site which could contribute directly to DNA binding in the NF-kappa-B complex.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286    E\_mail: info@targetmol.com    Address: 34 Washington Street, Wellesley Hills, MA 02481