

Anti-Acetyl-NF-kB p65 (Lys221) Polyclonal Antibody



Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 65 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity:

1. Western Blot analysis of HEPG2-UV cells using Acetyl-NFkB-p65 (K221) Polyclonal Antibody.
2. Western Blot analysis of HepG2-UV mouse brain cells using Acetyl-NFkB-p65 (K221) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized acetylpeptide: human NFκB-p65 around the acetylation site of K221

Antigen Species: Human

Synonyms: Acetyl-NF-κB p65 (K221);NF-κB p65K221 ac;Ac-NF-κB p65K221;Ac-NF-κB p65 (Lys221);NFκB3;p65;Ac-NF-κB p65 (K221);v-rel avian reticuloendotheliosis viral oncogene homolog A;NF-κB p65K221-acetylated

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481