


Anti-Acetyl-NFE2L2 (Lys599) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 68 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western Blot analysis of A549 cells using Acetyl-Nrf2 (K599) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized acetylpeptide: human Nrf2 around the acetylation site of K599
Antigen Species:	Human
Synonyms:	NFE2L2K599-acetylated;Ac-NFE2L2K599;NFE2L2K599ac;Acetyl-NFE2L2 (K599);Ac-NFE2L2 (K599);Ac-NFE2L2 (Lys599)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481