


## Anti-Phospho-NR4A1 (Ser351) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 64 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western Blot analysis of HepG2 cells using Phospho-Nur77 (S351) Polyclonal Antibody.



Application: ELISA,IHC,WB  
Recommended WB: 1:500-2000; IHC: 1:100-300;ELISA: 1:40000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: A synthesized phosphopeptide: human Nur77 around the phosphorylation site of S351  
Antigen Species: Human  
Uniprot ID: P22736  
Synonyms: NR4A1 (p-S351);NR4A1 (p-Ser351);p-NR4A1 (S351);p-NR4A1 (Ser351)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481