

Anti-p300 Antibody (7Z560)


Product Details

| | |
|-------------------|-----------------------|
| Reactivity: | Human,Mouse |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 300 kDa. |
| Clone: | 7Z560 |
| Purification: | Protein G |

Applications

Verified Activity:

1. Whole cell lysates of HeLa were separated by 10% SDS-PAGE, and the membrane was blotted with anti-p300 antibody.
2. Whole cell lysates of A549 were separated by 6% SDS-PAGE, and the membrane was blotted with anti-p300 antibody.



Application: ELISA,WB

Recommended WB: 1:500-2000; ELISA 1:5000-20000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide: human p300

Antigen Species: Human

Gene ID: 2033

Uniprot ID: Q09472

Synonyms: E1A-associated protein p300;Histone butyryltransferase p300;Histone acetyltransferase p300; Protein propionyltransferase p300;Protein lactyltransferas p300;p300 HAT;Histone crotonyltransferase p300;P300;Protein 2-hydroxyisobutyryltransferase p300;Ep300

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481