

## Anti-p300 Antibody (7Z560)


### Product Details

Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 300 kDa.
Clone:	7Z560
Purification:	Protein G

### Applications

#### Verified Activity:

1. Whole cell lysates of HeLa were separated by 10% SDS-PAGE, and the membrane was blotted with anti-p300 antibody.
2. Whole cell lysates of A549 were separated by 6% SDS-PAGE, and the membrane was blotted with anti-p300 antibody.



Application: ELISA,WB

Recommended WB: 1:500-2000; ELISA 1:5000-20000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

**Antigen Details**

Immunogen: A synthesized peptide: human p300  
Antigen Species: Human  
Gene ID: 2033  
Uniprot ID: Q09472  
Synonyms: E1A-associated protein p300;Histone butyryltransferase p300;Histone acetyltransferase p300; Protein propionyltransferase p300;Protein lactyltransferas p300;p300 HAT;Histone crotonyltransferase p300;P300;Protein 2-hydroxyisobutyryltransferase p300;Ep300

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481