

Anti-MAPK14 Antibody (1U459)


Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 42 kDa.
Clone:	1U459

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using P38 Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: P38
Uniprot ID:	Q16539
Synonyms:	Stress-activated protein kinase 2a;P38 MAPK;CSAID-binding protein;Mitogen-activated protein kinase 14;MXI2;EC 2.7.11.24;CSBP2;CSBP;MAPK 14;MAP kinase MXI2;Cytokine suppressive anti-inflammatory drug-binding protein;Mitogen-activated protein kinase p38 alpha;CSBP1;MAX-interacting protein 2;MAP kinase p38 alpha;MAP kinase 14;SAPK2A

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481