

Anti-MAPK3/ERK2 Antibody (7P251)



Product Details

| | |
|-------------------|--------------------------|
| Ig Type: | IgG1 |
| Reactivity: | Human,Rat,Mouse |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 44, 42 kDa. |
| Clone: | 7P251 |

Applications

Verified Activity:

1. Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain Tissue with P44/42 MAPK (ERK1/2) Mouse mAb diluted at 1:2,000.
2. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using P44/42 MAPK (ERK1/2) Mouse mAb diluted at 1:200.



| | |
|--------------|--------------------------------|
| Application: | IHC,WB |
| Recommended | WB: 1:1000-2000; IHC: 1:50-100 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

Immunogen: A synthesized peptide

Synonyms: MAPK 3;ERK 2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481