

## Anti-Acetyl-p53 (Lys305) Polyclonal Antibody


## Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 53 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

## Verified Activity:

1. Western Blot analysis of mouse-brain cells using primary antibody diluted at 1:1000(°C overnight). Secondary antibody Goat Anti-rabbit IgG IRDye 800( diluted at 1:5000, 25°C, 1 hour).



Application:	WB
Recommended	WB: 1:500-2000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized acetylpeptide: human p53 around the acetylation site of K305
Antigen Species:	Human
Synonyms:	Ac-p53 (K305);Acetyl-p53 (K305);Ac-p53 (Lys305);Ac-p53K305;p53K305-acetylated;p53K305ac

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481