

Anti-Acetyl-p53 (Lys320) Polyclonal Antibody



Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 53 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity:

1. Western Blot analysis of SKW3 cells using Acetyl-p53 (K320) Polyclonal Antibody.
2. Western Blot analysis of SKW3 cells using Acetyl-p53 (K320) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:20000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized acetylpeptide: human p53 around the acetylation site of K320
Antigen Species: Human
Synonyms: p53K320ac;p53K320-acetylated;Ac-p53 (K320);Ac-p53 (Lys320);Acetyl-p53 (K320);Ac-p53K320

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481