

## Anti-Acetyl-p53 (Lys372) Polyclonal Antibody



### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 53 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

### Applications

#### Verified Activity:

1. Western Blot analysis of 293T cells using Acetyl-p53 (K372) Polyclonal Antibody.
2. Western Blot analysis of Jurkat cells using Acetyl-p53 (K372) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:20000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

**Antigen Details**

Immunogen: A synthesized acetylpeptide: human p53 around the acetylation site of K372  
Antigen Species: human  
Synonyms: Ac-p53 (Lys372);p53K372-acetylated;Acetyl-p53 (K372);Ac-p53K372;Ac-p53 (K372);p53K372ac

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481