


Anti-Phospho-PAK1/2/3 (Ser144/141/139) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 65 kDa.
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis PAK1/2/3 (Phospho-Ser144+Ser141+Ser139) using HT-29 whole cell lysates.



Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human PAK1/2/3 around the phosphorylation site of Ser144+Ser141+Ser139
Antigen Species:	human
Uniprot ID:	Q13153 & Q13177 & O75914
Synonyms:	PAK1/PAK2/PAK3 (p-Ser144/Ser141/Ser139);p-PAK1/PAK2/PAK3 (Ser144/Ser141/Ser139); Phospho-PAK1/PAK2/PAK3 (S144/S141/S139);PAK1/PAK2/PAK3 (p-S144/S141/S139);p-PAK1/PAK2/PAK3 (S144/S141/S139)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481