


Anti-Phospho-PDK1 (Tyr243) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 49 kDa.
Purification:	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Applications

Verified Activity: 1. U87 cells were stimulated with or without hypoxia for 6 h. Immunoblotting analyses of mitochondrial lysates were performed with the indicated antibodies.



Application:	WB
Recommended	WB: 1:500-1000

Properties

Stability & Storage:	Store at 2°C-8°C for 1 month. Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Peptide sequence around phosphorylation site of tyrosine 243(D-L-Y(p)-Y-I) derived from Human PDHK1
Antigen Species:	human
Uniprot ID:	Q15118
Synonyms:	PDK1 (p-Tyr243);PDK1 (p-Y243);p-PDK1 (Tyr243);p-PDK1 (Y243)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481