


Anti-Phospho-PDLIM5 (Tyr251) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 64 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western Blot analysis of K562 cells using Phospho-PDLIM5 (Y251) Polyclonal Antibody.



Application: ELISA, WB
Recommended WB: 1:500-2000; ELISA: 1:10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human PDLIM5 around the phosphorylation site of Y251
Antigen Species: human
Uniprot ID: Q96HC4
Synonyms: PDLIM5 (p-Tyr251);p-PDLIM5 (Y251);p-PDLIM5 (Tyr251);PDLIM5 (p-Y251)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481