


Anti-Phospho-PERK (Ser1094) Polyclonal Antibody

Product Details

Ig Type: IgG
Reactivity: Human,Mouse,Rat
Conjugation: Unconjugated
Molecular Weight: Theoretical: 170 kDa.

Applications

Verified Activity: 1. Western blot analysis perk (Phospho-Ser1094) using Jurkat whole cell lysates.



Application: WB
Recommended WB: 1:1000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human perk around the phosphorylation site of Ser1094
Antigen Species: human
Uniprot ID: Q9NZJ5
Synonyms: p-PERK (Ser1094);p-PERK (S1094);PERK (p-Ser1094);PERK (p-S1094)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481