


Anti-Phospho-p53 (Thr55) Antibody (1F203)

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 53 kDa.
Clone:	1F203
Purification:	ProA affinity purified

Applications

Verified Activity: 1. Western blot analysis of Phospho-p53 (T55) on human skin lysates using anti-Phospho-p53 (T55) antibody at 1/1,000 dilution.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human p53 around the phosphorylation site of Thr55
Antigen Species:	human
Uniprot ID:	P04637
Synonyms:	p-p53 (T55);p-p53 (Thr55);p53 (p-Thr55);p53 (p-T55)

Research Background

p53, a DNA-binding, oligomerization domain- and transcription activation domain-containing tumor suppressor, upregulates growth arrest and apoptosis-related genes in response to stress signals, thereby influencing programmed cell death, cell differentiation, and cell cycle control mechanisms. p53 localizes to the nucleus, yet can be chaperoned to the cytoplasm by the negative regulator, MDM2. MDM2 is an E3 ubiquitin ligase that is upregulated in the presence of active p53, where it poly-ubiquitinates p53 for proteasome targeting. p53 fluctuates between latent and active DNA-binding conformations and is differentially activated through posttranslational modifications, including phosphorylation and acetylation. Mutations in the DNA-binding domain (DBD) of p53,

A DRUG SCREENING EXPERT

amino acids 110-286, can compromise energetically-favorable association with cis elements and are implicated in several human cancers.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481