

## Anti-Phospho-STAT3 (Tyr705) Antibody (2P792)


### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 79, 86 kDa.
Clone:	2P792

### Applications

#### Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Muscle Actin Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:100-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide
Uniprot ID:	P40763
Synonyms:	p-STAT3 (Tyr705);p-STAT3 (Y705);STAT3 (p-Y705);STAT3 (p-Tyr705)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481