

## Anti-Phospho-PI 3 Kinase Class 3 (Ser249) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 100 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

### Applications

Application:	WB
Recommended	WB: 1:1000-2000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human PI3 Kinase Class III around the phosphorylation site of Ser249
Antigen Species:	Human
Uniprot ID:	Q8NEB9
Synonyms:	p-PI 3 Kinase Class 3 (Ser249);p-PI 3 Kinase Class 3 (S249);PI 3 Kinase Class 3 (p-S249);PI 3 Kinase Class 3 (p-Ser249)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481