

Anti-Phospho-PRKCA/PKC beta (Thr638/641) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Actual: 76 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Applications

Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human PKC α around the phosphorylation site of Thr638 and 641
Antigen Species:	Human
Synonyms:	PRKCD;p-PRKCA (Thr638);PKC zeta;PRKCE;PKCG;p-PRKCA (T638);Protein kinase C alpha;PRKCA (p-T638);PKC epsilon;PKC alpha;PRKCB2;Protein kinase C gamma;PKC-A;Protein kinase C zeta;PKCA;Protein kinase C delta;PRK CZ;PKC2;Protein kinase C;PKCE;PRKCG;Phospho-PRKCA (T638);PRKCB1;PKCB;PKC beta;PKC gamma;PKC delta;PRKCA;Protein kinase C beta;PRKCA (p-Thr638);PRKCB;PKC-alpha;PKCD;KPCA_HUMAN;Protein kinase C epsilon

Research Background

protein kinase C alpha (PRKCA) Homo sapiens Protein kinase C (PKC) is a family of serine- and threonine-specific protein kinases that can be activated by calcium and the second messenger diacylglycerol. PKC family members phosphorylate a wide variety of protein targets and are known to be involved in diverse cellular signaling pathways. PKC family members also serve as major receptors for phorbol esters, a class of tumor promoters. Each member of the PKC family has a specific expression profile and is believed to play a distinct role in cells. The protein encoded by this gene is one of the PKC family members. This kinase has been reported to play roles in many different cellular processes, such as cell adhesion, cell transformation, cell cycle checkpoint, and cell volume control. Knockout studies in mice suggest that this kinase may be a fundamental regulator of cardiac contractility and Ca(2+) handling in myocytes. [provided by RefSeq, Jul 2]

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481