


Anti-Phospho-PKM (Tyr105) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 58 kDa.
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of extracts of various cell lines, using PKM2 (Phospho-Tyr105).



Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human PKM2 around the phosphorylation site of Tyr105
Antigen Species:	Human
Uniprot ID:	P14618
Synonyms:	CTHBP;OIP3;PKM2;pyruvate kinase, muscle;TCB;PKM2 (p-Y105);PKM2 (p-Tyr105);THBP1;p-PKM2 (Tyr105);HEL-S-30;PK3;p-PKM2 (Y105)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481