


## Anti-Plant actin Antibody (4N883)

### Product Details

Ig Type: IgG1  
Reactivity: Plant  
Conjugation: Unconjugated  
Molecular Weight: Theoretical: 43 kDa.  
Clone: 4N883

### Applications

Verified Activity: 1. Western blot analysis of Arabidopsis, diluted at 1) 1:5000 2) 1:10000.



Application: WB  
Recommended WB: 1:500-2000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Recombinant Protein

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481