


Anti-Plant actin Antibody (9X504)

Product Details

Ig Type:	IgG1
Reactivity:	Plant
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 43 kDa.
Clone:	9X504

Applications

Verified Activity: 1. Western blot analysis of Arabidopsis, (Q30) diluted at 1:5000.



Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481