

Anti-PR Antibody (1H838)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 38, 43 kDa.
Clone:	1H838
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Application:	IF,IHC
Recommended	IHC: 1:200; IF: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide
Uniprot ID:	P06401
Synonyms:	NR3C3;Progesterin receptor form A;Progesterin receptor form B;PRA;Nuclear receptor subfamily 3 group C member 3;PGR;PR;Progesterone receptor;PRGR;PRB

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481