

Anti-Phospho-PRC1 (Thr481) Polyclonal Antibody


Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 72 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Applications

Verified Activity:

1. Western blot analysis of 293T VEC lysate, antibody was diluted at 500. Secondary antibody was diluted at 1:20000.



Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Synthetic Phospho peptide from human protein at AA range: 481
Antigen Species:	Human
Uniprot ID:	O43663
Synonyms:	p-PRC1 (T481);p-PRC1 (Thr481);PRC1 (p-Thr481);PRC1 (p-T481)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481