


## Anti-Phospho-PKN1 (Thr774, 816) Polyclonal Antibody

### Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 104 kDa.
Purification:	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

### Applications

Verified Activity: 1. Western Blot analysis of lysates of 1. NIH/3T3 treated with LPS and 2. NIH/3T3, using primary antibody at 1:1000 dilution.



Application:	WB
Recommended	WB: 1:1000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized phosphopeptide: human PRK1/2 around the phosphorylation site of Thr774/816
Antigen Species:	Human
Uniprot ID:	Q16512
Synonyms:	PKN1 (p-T774, 816);p-PKN1 (Thr774, 816);PKN1 (p-Thr774, 816);p-PKN1 (T774, 816)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481